

Preserve with Confidence.  
Unleash Therapeutic Potential.

# CryoEase-PF

*Advanced MSC Cryopreservation  
Solution*

A photograph of two scientists in a laboratory setting. They are wearing white lab coats, blue surgical caps, and blue face masks. The scientist in the foreground is wearing white gloves and is using a pipette to transfer liquid from a multi-well plate. The scientist in the background is also wearing glasses and is looking at the plate. The background is slightly blurred, showing laboratory equipment and blue curtains.

CRYOEASE-PF

# ATLANTIS BIOSCIENCE PTE LTD

In Singapore, Atlantis Bioscience aims to enhance the lives of people worldwide through innovative life science solutions. We specialises in curating translational solutions for the research and development of cell & gene therapy, regenerative medicine, drug discovery & development, beauty & personal care, and food & nutritional science.

**‘Providing the finest **bench-to-bed support** to power you in making a difference.’**



*Engineered for Excellence,  
Trusted Globally*

# CryoEase-PF


hMSC DMSO & Protein Free Cryopreservation Media

## Overview

The **only DMSO-free, protein-free cryopreservation medium** that solves MSC storage challenges. Developed by **Singapore's cell therapy experts**, CryoEase-PF delivers superior viability and functionality while eliminating DMSO risks – setting the new standard for clinical-ready cell preservation.



## Key Features

- DMSO-free, protein-free formula
- Maintains 23-hour post thaw population doubling time (PDT)
- GMP-compliant materials 
- Preserve >85% Viability, Expansion Capacity & Potency after 72h at 4°C
- Dual-temperature stability



## How It Works

CryoEase-PF uses a proprietary, stress-free formula that:



Reduces ice  
crystal  
damage




Maintains  
osmotic  
stability



Shields cells from  
oxidative/metabolic  
stress

## A Fit for Every Frontline



### Clinical Settings

Supporting MSC-based  
therapies with  
predictable and  
reproducible outcomes.



### Biobanking

Enabling reliable  
experiments with high-  
quality stem cells.



### Research Laboratories

Enabling reliable  
experiments with high-  
quality stem cells.



### Veterinary Medicine

Facilitating innovative  
stem cell therapies for  
animal health.




# Real Results. Real Benefits.

**CryoEase-PF** isn't just a cryopreservation solution — it's a complete upgrade to how you store, transport, and use MSCs. Here's how it solves real pain points with real-world, tested results:


## ✓ Higher Cell Viability — Even After 72 Hours

MSCs stored in CryoEase-PF retained over 85% viability after 72 hours at 4°C, compared to less than 75% with leading alternatives.

**More viable cells = fewer repeats, lower costs, and stronger therapeutic or research outcomes.**



## ✓ Your Cells Stay Functional – Not Just "Alive"



MSCs stored in CryoEase-PF retain their ability to grow and expand, with doubling times of ~23 hours (freshly thawed) and ~30 hours (after 72 hours at 4°C).

**These are cells you can actually expand and work with — perfect for clinical therapies and reproducible lab experiments.**

✓ Preserves Stem Cell Identity and Potency

Flow cytometry histograms show that cells maintain expression of CD73, CD90, CD105 (positive markers), and low expression of CD11b, CD19, CD34, CD45, and HLA-DR (negative markers).





Table confirms 99%+ positivity for CD markers post-thaw and post-storage.





Microscopic images show robust adipogenic, osteogenic, and chondrogenic differentiation in both acutely thawed and 72hr stored MSCs.

Your cells stay true MSCs – reliable, reproducible, and compliant with ISCT standards.

# How **CryoEase-PF** Outperforms Competitors



Parameter	CryoEase-PF 	Commercial Alternatives 
<b>72h Viability</b>	>85%	<75%
<b>Population Doubling Time</b>	23h (acute thaw)	35–50h (variable)
<b>Toxicity Risks</b>	DMSO-free, xeno-free	DMSO-dependent
<b>ISCT Compliance</b>	Full retention	Often compromised



## Ready to Preserve Better?

Join researchers, clinicians, and biobanks worldwide who trust **CryoEase-PF** to protect the full potential of their stem cells.

Request a sample from us!  
[support@atlantisbioscience.com](mailto:support@atlantisbioscience.com)




Scan to order




# Step-by-Step Protocol

1



Trypsinise cells from monolayer.

2




Quantify cell number and viability.

3



Pellet cells by centrifugation.

4



Resuspend in an appropriate amount of **CryoEase-PF**.

5



Transfer to **cryovials**.

6



Equilibrate at room temperature for 10-15 mins.


*We recommend a freezing density of  $1-5 \times 10^6$  cells/mL*

7



Store at  $-80^{\circ}\text{C}$ .

*Optional: use control rate freezer/freezing devices like **CellHome** for  $-1^{\circ}\text{C}/\text{min}$  freezing*



# Discover leading edge solutions with us.



+65 8608 0974



[sales@atlantisbioscience.com](mailto:sales@atlantisbioscience.com)



[www.atlantisbioscience.com](http://www.atlantisbioscience.com)



visit our website

