# METHOD

# **Open Access**

# Automated digital image quantification of histological staining for the analysis of the trilineage differentiation potential of mesenchymal stem cells



Benjamin Eggerschwiler<sup>1,2</sup>, Daisy D. Canepa<sup>1,2</sup>, Hans-Christoph Pape<sup>1</sup>, Elisa A. Casanova<sup>1\*</sup> and Paolo Cinelli<sup>1,3\*</sup>

# Abstract

**Background:** Multipotent mesenchymal stem cells (MSCs) have the potential to repair and regenerate damaged tissues and are considered as attractive candidates for the development of cell-based regenerative therapies. Currently, there are more than 200 clinical trials involving the use of MSCs for a wide variety of indications. However, variations in their isolation, expansion, and particularly characterization have made the interpretation of study outcomes or the rigorous assessment of therapeutic efficacy difficult. An unbiased characterization of MSCs is of major importance and essential to guaranty that only the most suitable cells will be used. The development of standardized and reproducible assays to predict MSC potency is therefore mandatory. The currently used quantification methodologies for the determination of the trilineage potential of MSCs are usually based on absorbance measurements which are imprecise and prone to errors. We therefore aimed at developing a methodology first offering a standardized way to objectively quantify the trilineage potential of MSC preparations and second allowing to discriminate functional differences between clonally expanded cell populations.

**Method:** MSCs originating from several patients were differentiated into osteoblasts, adipocytes, and chondroblasts for 14, 17, and 21 days. Differentiated cells were then stained with the classical dyes: Alizarin Red S for osteoblasts, Oil Red O for adipocytes, and Alcian Blue 8GX for chondroblasts. Quantification of differentiation was then performed with our newly developed digital image analysis (DIA) tool followed by the classical absorbance measurement. The results from the two techniques were then compared.

**Result:** Quantification based on DIA allowed highly standardized and objective dye quantification with superior sensitivity compared to absorbance measurements. Furthermore, small differences between MSC lines in the differentiation potential were highlighted using DIA whereas no difference was detected using absorbance quantification.

**Conclusion:** Our approach represents a novel method that simplifies the laboratory procedures not only for the quantification of histological dyes and the degree of differentiation of MSCs, but also due to its color independence, it can be easily adapted for the quantification of a wide range of staining procedures in histology. The method is easily applicable since it is based on open source software and standard light microscopy.

**Keywords:** Digital image analysis, Dye quantification, Quantification of differentiation potential, Histology, Microscopy, Mesenchymal stem cells, Osteoblasts, Adipocytes, Chondroblasts

\* Correspondence: elisa.casanovazimmermann@usz.ch; paolo.cinelli@usz.ch

<sup>1</sup>Department of Trauma, University Hospital Zurich, Sternwartstrasse 14, 8091 Zurich, Switzerland

Full list of author information is available at the end of the article



© The Author(s). 2019 **Open Access** This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

# Background

Quantification of immunohistochemical staining by means of color analysis is a widely employed methodology in research and diagnostics. A major drawback of biochemical assays, such as absorbance measurement, is that it is destructive whereas quantification by image analysis is non-destructive. Furthermore, the possibility to perform automated analysis makes quantification faster, more objective, and less laborious than classical visual examination. Therefore, digital image analysis (DIA) has become a powerful and regularly used tool in a broad range of diagnostic and medical applications [1]. In pathology for example, DIA is used to analyze cell segmentation [2], mitosis [3], or nuclei [4]. Nevertheless, in other fields, this technology is not yet widely used. A possible field of application of DIA is the assessment of trilineage differentiation potential of mesenchymal stem cells (MSCs). Due to their low tumorigenic potential and their ability to rescue apoptotic cells after traumatic exposure, MSCs represent a promising source for regenerative medicine applications [5-8]. They were firstly identified in the bone marrow, but nowadays, they have been isolated and characterized from several adult and fetal tissues [9–12]. The most important challenge for the development of efficacious MSC-based therapies is currently the inability to consistently manufacture homogeneous populations of cells. MSC preparations do not consist of a pure population of cells but are rather a mixture of stem and linage-committed progenitor cells with heterogenous differentiation potential [13, 14]. The gold standard is the assessment of trilineage differentiation potential of MSCs to exhibit adipo-, chondro-, and osteogenesis as a measure of multipotency [15].

Various methods have been developed to characterize the degree of differentiation of MSCs based on transcriptome/proteome quantification. These readouts are very well established and quantitative. However, mRNA or protein levels do not necessarily reflect phenotypic features. The phenotypical characterization of such cell lines is often performed by histological staining [11, 12, 16–18] to prove differentiation down the adipogenic (Oil Red O staining), osteogenic (Alizarin Red S or von Kossa staining), and chondrogenic lineages (Alcian Blue 8GX or Safranin O staining). The main issue with histology is that the readout is binary. Histological staining confirms either the presence or absence of a specific cell type. Often only representative images are shown, lacking the confirmation that bright-field images are an accurate representation of the overall cell monolayer on the culture dish [16, 19-21]. In some cases, quantification of the staining is performed with absorbance measurements [22–25]. However, these protocols are laborious and vulnerable and show low sensitivity and inter-assay reproducibility. The fraction of stained cells is often vanishingly small compared to the total volume of solvent used to elute the dye from the cells. Furthermore, a considerable high amount of unspecifically bound dye is present leading to a low signal to noise ratio.

The development of quantifiable assays able to discriminate functional differences between populations and eventually predict the therapeutic efficacy of MSC products is urgently needed. Therefore, a method with a high signal to noise ratio and high sensitivity to quantify the degree of differentiation on a phenotypic level would be highly beneficial. We therefore aimed at developing a methodology offering first a standardized way to objectively quantify the trilineage potential of MSC preparations and second allowing to discriminate functional differences between clonally expanded cell populations. Moreover, this approach provides easily deployable metrics to assess differentiation efficacy toward osteoblasts, adipocytes, and chondroblasts using specific histological dyes (Alizarin Red S, Oil Red O, and Alcian Blue 8GX).

# Materials and methods

#### Gelatine coating

For gelatine coating,  $500 \ \mu l \ 0.1\% \ (w/\nu)$  gelatine in ddH<sub>2</sub>O (B. Braun, Germany, Cat. No. 0082479E) was added to each well. After 1 h incubation at room temperature, the gelatine solution was removed and plates were dried at room temperature.

# Cell isolation

Human MSCs were isolated from fat tissue with the consent of the patient according to the Swiss (KEK-ZH: StV 7-2009) and international ethical guidelines (ClinicalTrials.gov Identifier: NCT01218945) as reported previously [26]. The extraction procedure was performed like previously described [27]. MSCs were characterized according to established procedures. Of the 30 isolated primary MSCs, six were selected based on findings in a previous study concerning their differentiation capacity.

## Cell culture

Cells were cultured in DMEM (PAN Biotech, Germany, Cat. No. P04-03550) supplemented with 1% ( $\nu/\nu$ ) 100× Penicillin-Streptomycin Solution (Biowest, France, Cat. No. L0022), 10% ( $\nu/\nu$ ) FCS (Biowest, France, Cat. No. S181S), and 1% ( $\nu/\nu$ ) 200 mM L-glutamine solution (Sigma-Aldrich, USA, Cat. No. G7513), in a humidified atmosphere containing 5% CO<sub>2</sub>. The medium was exchanged every 3–4 days, and cells were split using 1× Trypsin-EDTA (Gibco/Life Technologies, USA, Cat. No. 25200-056) when reached 90–95% confluence.

For osteogenesis, MSCs were seeded on gelatinecoated Nunc<sup> $\infty$ </sup> 24-well plates (Thermo Fisher, USA, Cat. No. 142475) at a density of  $1.5 \times 10^4$  cells/cm<sup>2</sup>. After 24 h, differentiation was induced using StemPro<sup> $\circ$ </sup> Osteogenesis Kit (Gibco/Life Technologies, USA, Cat. No. A10072-01). Cells were regularly checked for morphology, and the medium was exchanged every 3–4 days.

For adipogenesis, MSCs were seeded on Nunc<sup>\*\*</sup> 48-well plates (Thermo Fisher, USA, Cat. No. 150687) at a density of  $2 \times 10^4$  cells/cm<sup>2</sup>. After 24 h, differentiation was induced using the StemPro<sup>\*</sup> Adipogenesis Kit (Gibco/Life Technologies, USA, Cat. No. A10070-01). Cells were regularly checked for morphology, and the medium was exchanged every 3–4 days.

For chondrogenesis, MSCs were seeded on Nunc<sup>\*\*</sup> 24-well plates (Thermo Fisher, USA, Cat. No. 142475) at a density of  $6 \times 10^3$  cells/cm<sup>2</sup>. After 24 h, differentiation was induced with the StemPro<sup>\*</sup> Chondrogenesis Kit (Gibco/Life Technologies, USA, Cat. No. A10071-01). Cells were regularly checked for morphology, and the medium was exchanged every 3–4 days.

## Cell staining

After 14, 17, and 21 days of differentiation, cells were fixed and stained for light microscopy.

Osteoblasts were washed with  $1 \times PBS$  (Kantonsapotheke Zürich, Switzerland, Cat. No. A171012) and fixed with 4% ( $\nu/\nu$ ) formaldehyde (Sigma, USA, Cat. No. F8775) in  $1 \times PBS$  for 30 min. After washing twice with ddH<sub>2</sub>O, Alizarin Red staining solution (0.7 g Alizarin Red S (Sigma, USA, Cat. No. A5533) diluted in 50 ml ddH<sub>2</sub>O at pH = 4.2) was added for 20 min. Afterwards, cells were washed four times with ddH<sub>2</sub>O, dried, and stored in the dark until image acquisition.

Adipocytes were washed with  $1 \times PBS$  and fixed with 10% ( $\nu/\nu$ ) formaldehyde in  $1 \times PBS$  for 65 min. After 5 min, the formaldehyde solution was exchanged. After washing twice with ddH<sub>2</sub>O, wells were rinsed with 60% ( $\nu/\nu$ ) 2-propanole (Sigma-Aldrich, USA, Cat. No. 59300) in ddH<sub>2</sub>O and dried. Oil Red O solution (0.15 g Oil Red O (Sigma-Aldrich, USA, Cat. No. 00625) diluted in 50 ml 60% ( $\nu/\nu$ ) 2-propanole in ddH<sub>2</sub>O) was added for 10 min. After four washing steps with ddH<sub>2</sub>O, the plates were dried and image acquisition was performed immediately after.

Chondroblasts were washed with  $1 \times PBS$  and fixed with  $4\% (\nu/\nu)$  formaldehyde in  $1 \times PBS$  for 20 min. Afterwards, cells were washed twice with ddH<sub>2</sub>O followed by a 3-min incubation at RT with  $3\% (\nu/\nu)$  acetic acid (Merck Millipore, Germany, Cat. No. 100063) in ddH<sub>2</sub>O. Alcian Blue solution (0.1 g Alcian Blue 8GX (Sigma, USA, Cat. No. A5268) diluted in 10 ml  $3\% (\nu/\nu)$  acetic acid in ddH<sub>2</sub>O at pH = 2.5) was added for 60 min followed by a washing step with 1 M HCl for 3 min. Afterwards, wells were washed four times with ddH<sub>2</sub>O. Multi-well dishes were afterwards stored in the dark until images were acquired.

# Absorbance measurement

Alizarin Red S was eluted from stained osteoblasts with 300 µl 10% (w/v) cetylpyridinium chloride in an aqueous 0.01 M Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub> solution at pH = 7 for 1 h. One hundred fifty microliters was transferred on a 96-well plate, and absorbance was measured at 560 nm. Ten percent (w/v) cetylpyridinium chloride in an aqueous 0.01 M Na<sub>2</sub>HPO<sub>4</sub>/NaH<sub>2</sub>PO<sub>4</sub> solution was used as blank.

Oil Red O was eluted from the stained adipocytes with 150  $\mu$ l 100% 2-propanole. After 1 min, the same amount of ddH<sub>2</sub>O was added and incubated for 4 min. One hundred fifty microliters was pipetted into a 96-well plate, and absorbance was measured at 520 nm. Fifty percent ( $\nu/\nu$ ) 2-propanole in ddH<sub>2</sub>O was used as blank.

### Image acquisition

Images were acquired with a Cytation 5 imaging reader (BioTek, USA, Cat. No. CYT5MPV). Whole wells were scanned in color bright-field mode with default settings regarding illumination and image capture. Single images were acquired with a horizontal overlap of 200  $\mu$ m and a vertical overlap of 250  $\mu$ m. One well on a 24-well plate needed 9 (horizontal) × 13 (vertical) images for stitching. For a 48-well plate, 7 (horizontal) × 10 (vertical) images were required.

#### Image processing and analysis

Single images were stitched together using Gen 5 image prime software (BioTek, USA, V3.03). Since the Cytation 5 imaging reader creates for each basic color (red, green, blue) a separate 16-bit image, a composite file was created using the open source software Cell Profiler [28]. Image analysis was performed with the open source software Fiji [29]. For any mathematical operation during the analysis steps, images were converted into a 32-bit floating point format to avoid pixel saturation. A macro was written in ImageJ macro language to analyze stitched images in batch mode. Analyzed parameters included area and area fraction.

#### Quantification of lineage-specific differentiation potential

For DIA, the measured values after 14, 17, and 21 days of differentiation were added up and then divided by 1000. For absorbance indeed, the measured values were added up and multiplied by 1000. This value is termed "Ranking Points".

#### Statistical analysis

All data analysis was carried out by using Graph Pad Prism 8 software (GraphPad Software, USA, Version 8.0.0). Unless otherwise stated, p values were calculated using Student's t tests and error bars indicate the standard deviation. Significance levels were subdivided into and \*\*\*\*(*p* value  $\leq$  0.0001). For correlation analysis, Pear-

\*(p value  $\leq 0.05$ ), \*\*(p value  $\leq 0.01$ ), \*\*\*(p value  $\leq 0.001$ ),

son's r value was calculated.

# Results

# Image acquisition and analysis protocol

In a first step, an image acquisition protocol was established to ensure consistent data generation. All relevant parameters, such as exposure time, white balance, light intensity, and plate dimensions, were fixed for all image acquisitions. The basic principle behind our method is the fact that a dye absorbs its complementary color leading to a specific bit depth value for each pixel (Fig. 1a). Our system is using a 16-bit CCD camera (2<sup>16</sup> bit depth increments), and samples are illuminated with three different LEDs with tight emission spectra and emission maxima at 460 nm (blue light), 523 nm (green light), and 623 nm (red light) to generate a color bright-field image (Fig. 1b). Because the relationship between spectral absorbance and stain concentration is only linear under monochromatic conditions, we can assume for our purposes that light transmission is linear with the pixel values in each channel [30, 31]. Therefore, the relative intensity of the histological dye is analyzed by dividing the signal of the dye by the signal of its complementary color. This means that for the red staining Alizarin Red S and Oil Red O, the red channel is divided by cyan, a mixture of green and blue, which represents the complementary color of red (Fig. 1c). A signal with a ratio of  $\frac{\text{red channel}}{\text{green channel+blue channel}} \ge 1$  was considered as a positive signal because the signal from the red channel is at least twice as strong as the average signal from the green and the blue channel. All pixels within the region of interest (ROI), which represents the bottom of the well with the cell monolayer, fulfilling this ratio were counted. Alizarin Red S and Oil Red O staining showed a very strong signal in the red channel. For Alizarin Red S, as an example, a mean pixel value of 34,437 was initially measured, indicating a low absorbance of the red light fraction. For the green channel, a mean pixel value of 4521 was determined and for the blue channel a mean pixel value of 5368. This indicates that most of the green and the blue light fraction is absorbed by the dye, leading to the red color perception. Because of the dominance of the red channel values, the green and the blue channel were neglected for the signal and the formula to define a positive pixel was maximally simplified. Oil Red O staining showed similar values. For Alcian Blue 8GX, which was used to confirm the presence of chondroblasts, the signal is found in the green and the blue channel and the corresponding complementary color, which is absorbed, is red. Therefore, a ratio of green channel+blue channel red -channel  $\geq 4$  represents a positive signal since the average signal intensity in the green and the blue channel is twice as high as in the signal in the red channel. Initial measurements of Alcian Blue 8GX stainings revealed that the green and the blue channel have similar mean pixel values (32,622 and 33,774) which are multiple times larger than the mean values from the red channel (9362). These values underline that mainly the red component of light is absorbed. Thus, the contribution of the red channel to the final signal was neglected, and the formula was again simplified as previously done for Alizarin Red S and Oil Red O quantification.

## Testing of the digital image analysis method

To test our approach, MSCs, isolated from the stromal vascular fraction of fat tissue from six healthy individuals that underwent adipose tissue excision, were differentiated into osteoblasts, adipocytes, and chondroblasts. After 14, 17, and 21 days of differentiation, cells were fixed and stained either with Alizarin Red S, Oil Red O, or Alcian Blue 8GX to confirm the corresponding cell lineage. Cells were cultured either in differentiation medium or in DMEM (control). Images of the entire well for each condition were taken for analysis (Fig. 1d). When we applied our quantification protocol (Additional file 1: Figure S1) to evaluate the differentiation in the MSC lines, we confirmed that the new DIA approach was able to reliably detect different color stains (Fig. 1e and Additional file 2: Figure S2). Moreover, independently from the dye used, the signal was nicely separated between differentiated and control cells and nearly all background signal was eliminated. (Fig. 1e). We were also able to detect changes during the differentiation into each of the three lineages and to precisely monitor the differentiation progress of MSCs over time (Fig. 1e).

Furthermore, with a × 40 magnification, our DIA approach allowed to analyze the histological staining on an extremely high resolution and the thresholding was highly effective. Noise was massively reduced and the whole system was resilient against external disturbing factors such as dirt. This was true for all the three staining types: Alizarin Red S, Oil Red O, and Alcian Blue 8GX (Fig. 1f).

# System validation—comparison with standard absorbance measurements

The current standard to quantify histological dyes is absorbance measurement. The dye needs to be eluted/extracted from the tissue, and then absorbance is measured at a dye-specific wavelength. In order to validate our system, we again differentiated MSCs into osteoblasts and adipocytes and stained the cells with the corresponding dye after 14, 17, and 21 days of differentiation. The same plate was first subjected to DIA



quantification, and thereafter the dye was extracted for absorbance measurements. Comparison of the two methods showed that DIA was more sensitive and eliminated most of the background signal (Fig. 2a). Cell line F14, as an example, that underwent osteogenic differentiation for 17 days showed a signal to noise ratio of nearly 19,000,000× with DIA, whereas with absorbance it was only  $48.6\times$  (Fig. 2a). Same applied for adipogenic differentiation. The background was almost completely eliminated with DIA leading to a signal to noise ratio of roughly 25,000× compared to a ratio of  $3.47\times$  from absorbance measurement (Fig. 2a).



represent SD). Six different cell lines were differentiated into osteogenic differentiation, n = 9 for adipogenic differentiation, end bass represent SD). Six different cell lines were differentiated into osteoblasts and adipocytes and afterwards analyzed with DIA and absorbance measurements (**b**). (n = 9 for all cell lines except for F22 n = 3, error bars represent SD). The ranking of all six cell lines showed a signal overestimation for absorbance measurements and a superior sensitivity of the DIA approach (**c**). Correlation analysis for osteogenesis revealed that DIA does reliably quantify Alizarin Red S. The correlation between absorbance measurements and DIA for adipocytes is less strong (**d**). The reason for this is a signal overestimation originating from absorbance measurements for weakly differentiated cells (e.g., controls). DIA also allows to quantify histological dyes that cannot be quantified by absorbance (**e**). Chondrogenic differentiation was quantified with DIA using Alcian Blue 8GX-stained cells (n = 9, error bars represent SD). Chondroblasts that were differentiated in 3D were cut into multiple cross sections, and positive pixels were compared to the total number of pixels of the spheroid (**f**) leading to a quantitative expression of the chondrogenic differentiation (**g**)

Next to the DIA protocol with increased sensitivity, we established a system that allows us to rank different mesenchymal stem cell lines according to their differentiation potential toward a distinct lineage. For therapeutic applications of MSCs, it is important to quantify the differentiation potential into the lineage of interest. For this purpose, we differentiated six MSC lines into osteoblasts and adipocytes (Fig. 2b). By adding up the measured values after 14, 17, and 21 days of differentiation, the differentiation potential for each cell line was expressed in a quantitative way and a classification of the lines was therefore possible (Fig. 2c). An ideal cell line produces a large amount of differentiated tissue in a short time period. Hence, our approach takes a time and a quantity component into consideration.

For osteogenic differentiation, it is clear that the line F10 and F27 have a poor ability to generate osteoblasts (Fig. 2b). Nevertheless, with our ranking approach, it is possible to observe that the line F27 differentiated better than the line F10 (Fig. 2c). These slight differences were not observed when absorbance was measured (Fig. 2c), thus confirming the superior sensitivity of the DIA approach compared to the classical absorbance approach. All other cell lines differentiated better and showed a similar differentiation pattern with DIA as well as with absorbance (Fig. 2c). Thus, the use of the highly sensitive DIA method in combination with our quantification procedure allows to highlight the slightest differentiation entern that offerentiating cell lines with very low signals.

Adipogenic differentiation was easily measured and quantified with DIA whereas absorbance measurements were extremely difficult, laborious, and vulnerable (Fig. 2b and Fig. 2c). Because Oil Red O staining is eluted with pure 2-propanole, a highly volatile solvent, measured absorbance values may vary massively. Already a slight difference in time between elution and absorbance measurement or variation in room temperature may lead to differences in evaporation and therefore differences in absorbance. Another drawback of Oil Red O absorbance measurement is the residual dye on the sidewall of the wells with cells that underwent differentiation. This residual dye is obviously also eluted leading to false results. Already the smallest amounts lead to unwanted signals in absorbance measurements and therefore compromise the detection limit and lead to a signal overestimation. Especially, cell lines that differentiate in a poor way, such as cell lines F4, F10, and F27, show heavily overestimated absorbance values since the true values are "covered" with the signal from the residual dye on the sidewalls (Fig. 2b and Additional file 3: Figure S3). Only for cell lines that differentiated well into adipocytes, the noise became negligible (Fig. 2b). Furthermore, absorbance measurements lead to wrong conclusions regarding the adipogenic differentiation potential (Fig. 2c). Thus, especially for adipogenic differentiation potential measurements, a highly sensitive method like DIA is of extreme importance, since the dye elution step, which is the one causing the measurement's variability, is completely omitted.

Correlation analysis between DIA and absorbance values revealed a very good correlation (Pearson's r = 0.9668) between DIA and absorbance measurement for Alizarin Red S quantification and a poor correlation (Pearson's r = 0.5438) for Oil Red O (Fig. 2d).

The weak correlation of Oil Red O is caused by the very limited detection capabilities of absorbance measurements which originate from the residual dye on the sidewalls. In our experiments, approximately 22,400 MSCs were exposed to adipogenic medium per well. Even if all cells differentiated into adipocytes, the total volume of Oil Red-stained lipids would still be extremely small compared to the total volume. The final volume during the elution process was 300  $\mu$ l. It is known that, on average, a single adipocyte has a volume of 0.0003  $\mu$ l, leading to a total volume of 6.72  $\mu$ l Oil Red O-stained lipids if all cells are differentiated into adipocytes [32]. These low volumes make absorbance measurements extremely susceptible to noise from unspecifically bound Oil Red O stain.

Osteogenic differentiation is less affected. Osteoblasts produce larger amounts of extracellular matrix which are stained by Alizarin Red S. Therefore, the residual dye on the sidewalls of each well is less disturbing.

Another advantage of the DIA method is its ability to quantify histological dyes that cannot be brought into a solution. Alcian blue is an example for such a dye. It is often used to show the presence of cartilage by staining the acidic carbohydrates which occur in the cartilage matrix. However, during standard staining protocols where Alcian Blue 8GX-stained samples are washed with water, Alcian Blue 8GX forms an insoluble phthalocyanine complex [33]. This is only soluble in highly concentrated sulfuric acid, which destroys the complex and therewith the chromatic properties of the dye. Therefore, no reliable protocol was established to measure Alcian Blue 8GX-stained cartilage from cell culture or tissue samples by absorbance spectroscopy. However, with our DIA approach, we were able to quantify the chondrogenic differentiation potential of six different MSC lines and found F22 to have the highest and F27 the lowest chondrogenic differentiation potential (Fig. 2e).

A major advantage of our DIA approach is its flexibility. It is color independent and not only applicable to cell monolayers in multi-well plates. As a proof of principle, chondroblasts grown in spheroids were also quantitatively assessed. Since all spheroids vary in their size, quantification via DIA was adjusted and the number of positive pixels was determined compared to the total number of pixels of the spheroid. Therefore, multiple cross sections of a spheroid that was grown in a chondrogenic medium were stained with Alcian Blue 8GX followed by image acquisition and analysis (Fig. 2f, g).

# Discussion

Over the past years with the increase of computational power, there has been significant progress in image analysis [34, 35]. Similarly, the mesenchymal stem cell research field expanded in the last years exponentially due to the unique therapeutic properties of these cells. Nevertheless, although the excitement on the potential use of MSCs in regenerative medicine is very big, the development of robust techniques for their unbiased characterization and classification is still very poor. We describe here the principle and the validation of an image analysis approach that simplifies the laboratory procedures to objectively quantify and classify the degree of differentiation as well as the differentiation potential among different MSC cell lines or cell subpopulations. This approach is of course not limited to MSCs only but can be applied to any type of cells and tissue and also staining dye.

The method proposed here produces robust results with a superior sensitivity compared to absorbance measurements. In fact, very low signals, such as Oil Red O-stained poorly differentiated adipocytes, could be detected and analyzed with our approach whereas absorbance measurement was only reliable for large values of very well-differentiated cell lines.

Furthermore, absorbance measurements are prone to errors, and the sample preparation procedures are time-consuming and difficult to standardize. Since Oil Red O is eluted with pure 2-propanole, a highly volatile solvent, a lot of variability is added to the absorbance measurement and makes the whole measurement highly vulnerable to noise. All Oil Red O and Alizarin Red S stainings of differentiated cells lead to residual dye on the wall of a well, lowering the sensitivity and signal to noise ratio when quantified by absorbance measurements. With our approach, we completely eliminate the distaining step, avoiding the introduction of additional noise to the measurements.

On the other hand, variance is increased in DIA compared to absorbance measurements. Due to the high sensitivity, little differences which are not detected by absorbance measurements may lead to a considerable signal difference during DIA quantification.

Compared to other image analysis procedures, such as color deconvolution, our system has a superior sensitivity. Color deconvolution is based on 8-bit RGB images and needs a preliminary definition of the color deconvolution vectors. This definition is difficult and therefore prone to errors. Our bit-depth-independent DIA approach works with 16-bit images which represents a  $256 \times$  fold increase in sensitivity/resolution compared to 8-bit images, and the difficult vector definition step is omitted. Moreover, color deconvolution is limited to stoichiometric staining reactions [36]. Thus, color deconvolution is not suitable to

quantify immunohistochemical stainings. Our approach also allows to quantify such stainings because we only consider pixels above or below a predefined threshold either as positive or negative. We do not quantify the spatial amount of staining as it is done during color deconvolution.

# Conclusion

The new proposed method clearly proved its superiority compared to the conventionally used methods for MSC differentiation, quantification, and classification, and it represents therefore for the MSC research field a very important new tool.

Moreover, a major advantage of this method is its flexibility: Neither it is limited to a certain dye or cell type nor to cell monolayers in multi-well plates, and by using *Fiji's* batch mode, it is possible to automatically analyze large numbers of sections. Data credibility is also increased with availability of the original images.

In conclusion, with its high sensitivity, our approach can be used as a diagnostic tool because it allows to quantitatively distinguish cell lines and their differentiation capability based on their phenotype.

# **Additional files**

Additional file 1: Figure S1. Workflow for DIA. The flow chart shows all crucial steps for proper image analysis/quantification. Each step is described in the "Materials and methods" section in detail. (PDF 296 kb)

Additional file 2: Figure S2. Illustration of original and thresholded images. In this graph, original images are shown together with its corresponding thresholded image for all three histological dyes used, Alizarin Red S for osteoblasts, Oil Red O for adipocytes, and Alcian Blue 8GX for chondroblasts. The yellow circle highlights the ROI in which the pixels were analyzed. (PDF 493 kb)

Additional file 3: Figure S3. Illustration of superiority of the DIA approach. The images illustrate the differentiation into adipocytes after 14, 17, and 21 days of differentiation. Cell line F14 represents a cell line with a high adipogenic potential whereas F27 represents a cell line with a low adipogenic potential. The comparison between DIA and absorbance measurements reveals that absorbance measurements overestimate low signals. According to absorbance measurements, the degree of differentiation of cell line F27 after 21 days is roughly the same as for cell line F14 after 17 days of differentiation. Absorbance measurements are extremely prone to errors because a lot of unspecifically bound dye (especially from the side walls of a well) is brought into the solution. (PDF 488 kb)

#### Abbreviations

CCD: Charge coupled device; DIA: Digital image analysis; MSC: Mesenchymal stem cells; RGB: 3 channels (red, green, blue); ROI: Region of interest

#### Acknowledgements

We thank Sonja Hemmi and Sonja Märsmann for their assistance in the laboratory during the entire project. Image accquisition and analysis was performed with equipment and support of the Center for Microscopy and Image Analysis, University of Zurich. Further we thank the Epidemiology, Biostatistics and Prevention Institute for their statistical consulting.

#### Funding

This research was funded by the Department of Trauma, University Hospital Zurich.

#### Availability of data and materials

Please contact the authors for data requests.

#### Authors' contributions

BE performed all experiments and data analysis, contributed to the experimental design, and wrote the computer scripts and the manuscript. DDC performed the experiments. EAC, H-CP, and PC directed the study and wrote the manuscript. All authors read and approved the final manuscript.

#### Ethics approval and consent to participate

Human cells were isolated with the consent of the patient according to Swiss (KEK-ZH: StV 7-2009) and international ethical guidelines (ClinicalTrials.gov Identifier: NCT01218945).

#### Consent for publication

Not applicable.

#### Competing interests

The authors declare that they have no competing interests.

#### **Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

#### Author details

<sup>1</sup>Department of Trauma, University Hospital Zurich, Sternwartstrasse 14, 8091 Zurich, Switzerland. <sup>2</sup>Life Science Zurich Graduate School, University of Zurich, Winterthurerstrasse 190, 8057 Zurich, Switzerland. <sup>3</sup>Center for Applied Biotechnology and Molecular Medicine, University of Zurich, Winterthurerstrasse 190, 8057 Zurich, Switzerland.

# Received: 7 December 2018 Revised: 9 January 2019 Accepted: 11 February 2019 Published online: 26 February 2019

#### References

- Prasad K, Prabhu GK. Image analysis tools for evaluation of microscopic views of immunohistochemically stained specimen in medical research-a review. J Med Syst. 2012;36(4):2621–31.
- Song YY, Tan EL, Jiang XD, Cheng JZ, Ni D, Chen SP, et al. Accurate cervical cell segmentation from overlapping clumps in pap smear images. leee T Med Imaging. 2017;36(1):288–300.
- Albarqouni S, Baur C, Achilles F, Belagiannis V, Demirci S, Navab N. AggNet: deep learning from crowds for mitosis detection in breast cancer histology images. leee T Med Imaging. 2016;35(5):1313–21.
- Kashif MN, Raza SEA, Sirinukunwattana K, Arif M, Rajpoot N. Handcrafted features with convolutional neural networks for detection of tumor cells in histology images. In Proc. IEEE Int. Symp. Biomed. Imag., 2016;1029–32.
- Bhartiya D. Are mesenchymal cells indeed pluripotent stem cells or just stromal cells? OCT-4 and VSELs biology has led to better understanding. Stem Cells Int. 2013;2013:547501.
- Cselenyak A, Pankotai E, Horvath EM, Kiss L, Lacza Z. Mesenchymal stem cells rescue cardiomyoblasts from cell death in an in vitro ischemia model via direct cell-to-cell connections. BMC Cell Biol. 2010;11.
- Kim SY, Lee JH, Kim HJ, Park MK, Huh JW, Ro JY, et al. Mesenchymal stem cell-conditioned media recovers lung fibroblasts from cigarette smokeinduced damage. Am J Physiol-Lung C. 2012;302(9):L891–908.
- Li N, Sarojini H, An J, Wang E. Prosaposin in the secretome of marrow stroma-derived neural progenitor cells protects neural cells from apoptotic death. J Neurochem. 2010;112(6):1527–38.
- Ayatollahi M, Geramizadeh B, Zakerinia M. Human bone marrow-derived mesenchymal stem cell: a source for cell-based transplantation. Transpl Int. 2013;26:185–6.
- Ciavarella S, Stucci S, De Matteo M, Silvestris F. Umbilical cord stroma as source of mesenchymal stem cells for bone regeneration. Bone. 2009;45: S143–S4.

- In 'tAnker PS, Scherjon SA, Kleijburg-van der Keur C, Noort WA, FHJ C, Willemze R, et al. Amniotic fluid as a novel source of mesenchymal stem cells for therapeutic transplantation. Blood. 2003;102(4):1548–9.
- Russo V, Yu C, Belliveau P, Hamilton A, Flynn LE. Comparison of human adipose-derived stem cells isolated from subcutaneous, omental, and intrathoracic adipose tissue depots for regenerative applications. Stem Cell Transl Med. 2014;3(2):206–17.
- Mo M, Wang S, Zhou Y, Li H, Wu Y. Mesenchymal stem cell subpopulations: phenotype, property and therapeutic potential. Cell Mol Life Sci. 2016; 73(17):3311–21.
- Pittenger MF, Mackay AM, Beck SC, Jaiswal RK, Douglas R, Mosca JD, et al. Multilineage potential of adult human mesenchymal stem cells. Science. 1999;284(5411):143–7.
- Dominici M, Le Blanc K, Mueller I, Slaper-Cortenbach I, Marini FC, Krause DS, et al. Minimal criteria for defining multipotent mesenchymal stromal cells. The International Society for Cellular Therapy position statement. Cytotherapy. 2006;8(4):315–7.
- Gharibi B, Hughes FJ. Effects of medium supplements on proliferation, differentiation potential, and in vitro expansion of mesenchymal stem cells. Stem Cell Transl Med. 2012;1(11):771–82.
- Stanko P, Kaiserova K, Altanerova V, Altaner C. Comparison of human mesenchymal stem cells derived from dental pulp, bone marrow, adipose tissue, and umbilical cord tissue by gene expression. Biomed Pap. 2014; 158(3):373–7.
- Taghiyar L, Hosseini S, Hesaraki M, Sayahpour FA, Aghdami N, Eslaminejad MB. Isolation, characterization and osteogenic potential of mouse digit tip Blastema cells in comparison with bone marrow-derived mesenchymal stem cells in vitro. Cell J. 2018;19(4):585–98.
- Chan CKF, Gulati GS, Sinha R, Tompkins JV, Lopez M, Carter AC, et al. Identification of the human skeletal stem cell. Cell. 2018;175(1):43–56 e21.
- Desance M, Contentin R, Bertoni L, Gomez-Leduc T, Branly T, Jacquet S, et al. Chondrogenic differentiation of defined equine mesenchymal stem cells derived from umbilical cord blood for use in cartilage repair therapy. Int. J. Mol. Sci. 2018;19:537.
- Ikebuchi Y, Aoki S, Honma M, Hayashi M, Sugamori Y, Khan M, et al. Coupling of bone resorption and formation by RANKL reverse signalling. Nature. 2018;561(7722):195.
- 22. Bjornsson S. Quantitation of proteoglycans as glycosaminoglycans in biological fluids using an alcian blue dot blot analysis. Anal Biochem. 1998;256(2):229–37.
- Frazier SB, Roodhouse KA, Hourcade DE, Zhang L. The quantification of glycosaminoglycans: a comparison of HPLC, carbazole, and Alcian Blue methods. Open Glycosci. 2008;1:31–9.
- Newton PT, Staines KA, Spevak L, Boskey AL, Teixeira CC, Macrae VE, et al. Chondrogenic ATDC5 cells: an optimised model for rapid and physiological matrix mineralisation. Int J Mol Med. 2012;30(5):1187–93.
- Yu F, Dong Y, Yang YW, Lin PT, Yu HH, Sun X, et al. Effect of an experimental direct pulp-capping material on the properties and osteogenic differentiation of human dental pulp stem cells. Sci Rep-Uk. 2016;6.
- Buschmann J, Gao SP, Harter L, Hemmi S, Welti M, Werner CML, et al. Yield and proliferation rate of adipose-derived stromal cells as a function of age, body mass index and harvest site-increasing the yield by use of adherent and supernatant fractions? Cytotherapy. 2013;15(9):1098–105.
- Zuk PA, Zhu M, Mizuno H, Huang J, Futrell JW, Katz AJ, et al. Multilineage cells from human adipose tissue: implications for cell-based therapies. Tissue Eng. 2001;7(2):211–28.
- Carpenter AE, Jones TR, Lamprecht MR, Clarke C, Kang IH, Friman O, et al. CellProfiler: image analysis software for identifying and quantifying cell phenotypes. Genome Biol. 2006;7(10).
- Schindelin J, Arganda-Carreras I, Frise E, Kaynig V, Longair M, Pietzsch T, et al. Fiji: an open-source platform for biological-image analysis. Nat Methods. 2012;9(7):676–82.
- Haub P, Meckel T. A model based survey of colour deconvolution in diagnostic brightfield microscopy: error estimation and spectral consideration. Sci Rep. 2015;5:12096.
- Ruifrok AC, Johnston DA. Quantification of histochemical staining by color deconvolution. Anal Quant Cytol Histol. 2001;23(4):291–9.
- Leonhardt W, Hanefeld M, Haller H. The adipocyte volume in human adipose tissue: 1. Lipid space, normal and maximum values, and the relation to body weight index. Int J Obes. 1978;2(1):33–45.
- Kiernan JA. Dyes and other colorants in microtechnique and biomedical research. Color Technol. 2006;122(1):1–21.

- Eklund A, Dufort P, Forsberg D, LaConte SM. Medical image processing on the GPU - past, present and future. Med Image Anal. 2013;17(8):1073–94.
- Riber-Hansen R, Vainer B, Steiniche T. Digital image analysis: a review of reproducibility, stability and basic requirements for optimal results. Apmis. 2012;120(4):276–89.
- van der Loos CM. Multiple immunoenzyme staining: methods and visualizations for the observation with spectral imaging. J Histochem Cytochem. 2008;56(4):313–28.

# Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

#### At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

